

RECIRCULATOR

ART. NO. 150010

GENERAL

The KEOFITT RECIRCULATOR is a device to enable chemical cleaning/disinfection of KEOFITT sampling valves and KEOFITT sampling bags.

Compatible to all KEOFITT sampling valves with W9, 4KAI, Basix and Sesame hose piece as well as the KEOFITT sampling bags.

Cleaning agents not included.

CAUTION! As with any chemical solutions the care, handling and compatibility with the equipment must be considered before undertaking use of the device for sample valve cleaning. Follow safety measures to avoid chemical contact.



FEATURES

Connection: W9 in-/outlet and KEOFITT Sampling Bag outlet

Inlet / Outlet: Quick Coupling

CERTIFICATION*

O-ring: EU EC 1935/2004 · FDA CFR 21 §177.2600 · USP Class VI · ADI Free

TECHNICAL DATA

Material (process contact)

-

Material (without process contact)



Bottle (Capacity: 1000 ml)	Polypropylene (PP)
Screw cap	Polypropylene (PP)
Sealing ring	Thermoplastic elastomer (TPE)
Hoses	PVC
Check valves	Nylon (PA)
Quick Coupling	Stainless Steel AISI 316L (1.4404)
O-ring	Silicone, Red (#800830)

Net weight

Kg/lbs: 0.465 kg / 1.024 lbs

SPARE PARTS

#800830	10 pack O-ring EPDM 10.3 x 2.4
#150089	QC KSB/W9 hose barb Ø4
#550311	Bulb pump
#550314	steel top
#800088	QC W9 hose barb Ø4



*For further information and downloads please visit www.keofitt.dk