



## CLIP-ON STEAMER ART. NO. 550009

## **GENERAL**



The KEOFITT CLIP-ON STEAMER is a steam source that offers the ability to steam your KEOFITT sampling valve where steam is otherwise unavailable. All you need is water and a torch.



To be used with all steam sterilizable sampling valves.



Torch is NOT included!



CAUTION! Open flame, it is recommended to use gloves during the sterilising process as both CLIP-ON STEAMER, fittings and sampling valve become hot during the process

# **TECHNICAL DATA**

### **Material (without product contact)**

· Steel parts Stainless steel

· O-ring Silicone

#### **Capacity**

· Capacity 50 ml


**Net Weight** 

· Kg/lbs 0,500 kg

## **Spareparts**

800830 10 PACK 0-RING 10,3 x 2,4 SILICONE FDA





Data Sheet from keofitt.dk

Last updated 09-02-2017