







RECIRCULATOR

ART. NO. 550010

GENERAL

-  The KEOFITT RECIRCULATOR is a device to enable chemical cleaning/disinfection of KEOFITT sampling valves.
-  Compatible to all KEOFITT sampling valves with W9, 4KAI, Basix and Sesame hose piece.
-  Cleaning agents not included.
-  **CAUTION!** As with any chemical solutions the care, handling and compatibility with the equipment must be considered before undertaking use of the device for sample valve cleaning. Follow safety measures to avoid chemical contact.

FEATURES

-  Connection: W9 hose piece in-/outlet
-  Inlet / Outlet: Quick Coupling

TECHNICAL DATA

Material (process contact)

Material (without process contact)

- Bottle Polypropylene (PP)
- (Capacity: 1000 ml)
- Screw cap Polypropylene (PP)
- Sealing ring Thermoplastic elastomer (TPE)
- Hoses PVC
- Check valves Nylon (PA)
- Quick Coupling Stainless Steel AISI 316L (1.4404)
- O-ring Silicone, RED (#800830)
- EU EC 1935/2004 · FDA CFR 21 §177.2600
- USP Class VI · ADI Free

Net Weight

- kg/lbs 0.465 kg/1.024 lbs

Spareparts

- 800830 10 PACK O-RING EPDM 10.3 x 2.4

